FleXNeedle®

REF Catalog number 10005 (18 gauge)

Instructions For Use (IFU)


Figure 1: FlexNeedle with the needle retracted and stylet removed.

WARNINGS AND PRECAUTIONS

4 Warnings / Precautions

- The FlexNeedle should only be used by physicians thoroughly trained in bronchoscopy and in the particular technique and procedure to be performed, and familiar with the associated risks. These operating instructions must be read, understood, and followed.
- The FlexNeedle should be used with caution and only after careful consideration on patients who are at risk of bleeding complications.
- If manually thrusting the catheter to penetrate, **DO NOT USE** excessive force as this may create a passage of uncontrolled depth that could result in serious harm to the patient, such as excessive bleeding or pneumothorax.
- The FlexNeedle should only be used in a manner consistent with the Bronchoscope manufacturer's Instructions for Use. Use of the FlexNeedle in contradiction to the instructions may harm the device, bronchoscope or patient.

5 Possible Complications

- Bleeding
- Pneumomediastinum
- Pneumothorax
- Infection

6 Preparation for Use

6.1 Materials Required

- Broncus FlexNeedle
- Luer-lock aspiration syringe
- Flexible bronchoscope with a working channel of 2.0 mm or greater

6.2 Device Inspection and Preparation

- 6.2.1 Read all IFUs, package inserts, labels, and warnings for the FlexNeedle before beginning the procedure.

CAUTION: Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.

Table of Contents

- 1 Device Description
- 2 Intended Use
- 3 Contraindications
- 4 Warnings / Precautions
- 5 Possible Complications
- 6 Preparation for Use
- 7 Operating Instructions
- 8 Storage

1 Device Description

The FleXNeedle is a catheter with an 18 gauge coring needle tip in a protective sheath that fits down the 2-mm working channel of a standard flexible optical bronchoscope. It is intended for transbronchial retrieval of tissue samples.

The FleXNeedle is comprised of a handle, needle, sheath and removable stylet. The handle actuates the needle and allows for aspiration and the optional use of a stylet. The sheath protects the bronchoscope and user from puncture by the needle. The FleXNeedle is supplied sterile, for single-patient-use.

Product Configuration and Dimensions

The FleXNeedle (see illustration in Figure 1) has the following specifications:

Catalog Number	Needle Penetration Depth	Maximum Catheter Outer Diameter (OD)	Needle Penetration Depth
10005	18 ga	1.9 mm	0 - 20 mm, in 5 mm increments

2 Intended Use

The FleXNeedle is indicated for use in aspiration in carinal, paratracheal and hilar lesions of the bronchial tree where biopsy forceps cannot obtain a submucosal sample.

3 Contraindications

There are no known contraindications.

- 6.2.2 Aseptically remove the FleXNeedle from the package and inspect for any damage, such as broken or crushed areas of the catheter shaft, sharp or protruding edges at the distal tip, or any kinks in the catheter shaft.

CAUTION: If any damage or irregularity is found, **DO NOT USE** the device as this could result in harm to the patient.

CAUTION: DO NOT USE the device past its "Use Before" (expiration) date.

- 6.2.3 Before inserting the device into the working channel of the bronchoscope, visually confirm that the needle is retracted and locked in the "needle-in" position within the sheath.

CAUTION: The needle component of the FleXNeedle should be retracted within the sheath at all times, except during airway wall penetration and tissue aspiration. Inserting a device with an exposed needle into the working channel may cause damage to the bronchoscope and could result in serious harm to the patient.

CAUTION: DO NOT place or manipulate the catheter in the airways unless under direct visualization with a bronchoscope as this may result in harm to the patient.

7 Operating Instructions

- 7.1 Confirm the needle is locked in the "needle-in" position (see Figure 1). If using the stylet, confirm the stylet luer cap is firmly tightened on the luer fitting of the handle.
- 7.2 Carefully advance the FleXNeedle through the working channel of the bronchoscope using short (approximately 20 mm) increments, until the distal end of the catheter is within bronchoscopic view.

CAUTION: Always insert and remove this device with the articulating portion of the bronchoscope in the neutral position to avoid damage to the device and/or bronchoscope.

CAUTION: Use care when handling the device and **DO NOT** use excessive force while advancing or manipulating the catheter in the bronchoscope. This can result in (a) kinks and/or leaks in the aspiration lumen of the needle, (b) damage to the bronchoscope, or (c) harm to the patient's airway wall.

CAUTION: If any damage occurs, including a kink or an aspiration leak, **DO NOT USE** and discard.

- 7.3 Advance the catheter to a position distal of the bronchoscope. Contact the airway wall with the sheath and maintain position in preparation for penetration.
- 7.4 Place the needle actuation stop at the desired depth. Penetrate the airway wall by positioning the sheath against airway wall. Then use the handle thumb ring to firmly thrust the needle tip to puncture the airway wall.
- CAUTION:** Use care to avoid cartilage rings when attempting to penetrate airway walls, to prevent needle penetration failure and bending. Bending may damage the needle.
- 7.5 For sample collection, if using a stylet, remove it. Attach a luer-lock aspiration syringe to the luer fitting on the FleXNeedle. Place the needle actuation stop to define your desired needle penetration depth (up to 20 mm).

CAUTION: DO NOT apply advancing pressure on the thumb ring while moving the actuation stop. This can result in the needle advancing to a depth that is further than intended.

- 7.6 Use the handle thumb ring to firmly thrust the needle tip to the desired depth, to capture and enclose a sample. Hold the needle extended using the thumb ring, and apply suction with the aspiration syringe.

7.6.1 To oscillate the needle in the target, use the handle and the needle penetration depth positions to retract and advance the needle the desired amount.

- 7.7 Withdraw the needle from the target tissue and release the vacuum.

- 7.8 Retract the needle by withdrawing the thumb ring on the handle until the needle is completely within the sheath. If the needle does not completely retract initially, apply a rapid withdrawal action using the handle thumb ring. Additionally, if necessary, relax the bronchoscope articulation and repeat the rapid withdrawal action using the handle thumb ring. Visualize the device tip to confirm that the needle has been fully retracted; lock the needle in the "needle-in" position. Remove the FleXNeedle from the bronchoscope.

CAUTION: Removing the device without retracting the needle completely into the sheath may result in damage to the bronchoscope.

- 7.9 Expel the sample with air or saline, or push with the stylet, as desired. Proceed to handle the sample per your institution's preferred method.
- 7.10 For additional passes, verify device integrity, reinsert the stylet if desired, and repeat above steps.
- 7.11 At the completion of the patient procedure, discard the FleXNeedle in the appropriate disposal container.

8 Storage

Store at controlled room temperature.

LEGAL NOTICE

© 2013 Broncus Medical, Inc. All rights reserved. This product is covered by one or more of the following patents: U.S. Patent No. 8,517,955. Patents Pending. BRONCUS, the BRONCUS logo, and FleXNeedle are registered trademarks of Broncus Medical, Inc.

Graphic Symbol Legend for Medical Device Labeling

	Catalog number
	Federal law (USA) restricts this device to sale by or on the order of a physician.
	Consult instructions for use
	Irradiation sterilized. Sterility guaranteed if package unopened and undamaged.
	For single-patient-use only.
	Use by
	For single-patient-use only.
	Do not use if package is damaged
	Manufacturer
	European Representative

	Manufacturer: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Suite A8 Mountain View, CA, 94043 USA (877) 428-1600
	European Representative: C.K. Gill, Medical Consulting Plus PO Box 3112 6202 NC Maastricht, Netherlands +31-43-325-6771

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Produttore: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Suite A8 Mountain View, CA, 94043 USA (877) 428-1600