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1 About these Instructions for Use

1.1 Symbols Glossary

Symbol	Description
	Caution: Consult Instructions for Use
	Do not use if package is damaged
	Keep away from direct sunlight
	Keep dry
REF	Catalog number
LOT	Batch code
QTY	Quantity per packaging unit
	Manufacturer

Table 1: Symbols Glossary

1.2 Abbreviations

- WD: Washer Disinfectant

1.3 Terminology

- Distal: From the surgeon's view further away
- Proximal: From the surgeon's view closer

1.4 Safety Information Marking

WARNING

Non-observance may result in serious injuries, serious deterioration of the general condition or the death of the patient, user, or a third party.

NOTICE

Product damage or other damage may occur in case of non-compliance.

2 Important Safety Information

WARNING

- Before using the product, read the Instructions for Use. Adhere to and save the Instructions for Use.
Otherwise there are risks to the health of your patient.
- Application only by a physician trained in the procedure.
Otherwise there are risks to the health of your patient.
- Do not modify the product.
Otherwise there are risks to the health of your patient.

NOTICE

- Handle the product with utmost care. This applies especially to blades, fine tips, and other delicate areas.
Otherwise, the product could be damaged.

3 Product Codes / REF

REF	Name	For ROTEPS SHAFT
02BRT10200	Handle without ratchet	Ø 2 mm
02BRT10202	Handle with shiftable ratchet	Ø 2 mm
02BRT10400	Handle without ratchet	Ø 4 mm
02BRT10401	Handle with ratchet	Ø 4 mm
02BRT10402	Handle with shiftable ratchet	Ø 4 mm

Table 2: ROTEPS HANDLE

REF	A [cm]	B [mm]
02BRT20260	60	2
02BRT20460	60	4

Table 3: ROTEPS SHAFT

REF	Name	Jaw, closed Ø [mm]	For ROTEPS SHAFT
02BRT30294	Retraction forceps, atraumatic inner teeth, one movable jaw	3	Ø 2 mm
02BRT30295	Fenestrated forceps	3	Ø 2 mm
02BRT30296	Hook scissors	3	Ø 2 mm
02BRT30299	Grasping forceps, extra strong	3	Ø 2 mm
02BRT30497	Retraction forceps, atraumatic inner teeth	5	Ø 4 mm
02BRT30498	Y-Stent forceps	5	Ø 4 mm

REF	Name	Jaw, closed Ø [mm]	For ROTEPS SHAFT
02BRT30499	Grasping forceps, extra strong	5	Ø 4 mm

Table 4: ROTEPS HEAD

REF	Name	Comment
02BRT40010	Cleaning adapter for ROTEPS [®] SHAFT, with Luer Lock adapter incl. O-Ring	Accessories
02BRT40040	O-Ring, 0.7 x 1 mm, for cleaning adapter	Spare part

Table 5: ROTEPS CLEANSE

4 Scope of Delivery

Product	Package contents
ROTEPS HANDLE	<ul style="list-style-type: none"> • 1 x handle • 1 x Instructions for Use
ROTEPS SHAFT	<ul style="list-style-type: none"> • 1 x shaft • 1 x Instructions for Use
ROTEPS HEAD	<ul style="list-style-type: none"> • 1 x instrument attachment • 1 x Instructions for Use
ROTEPS CLEANSE (rinsing adapter)	<ul style="list-style-type: none"> • 1 x rinsing attachment • 1 x Instructions for Use
ROTEPS CLEANSE (O-Ring)	<ul style="list-style-type: none"> • 10 x O-Ring

Table 6: Package contents

5 Intended Use

For manipulating, grasping and cutting anatomic structures and foreign bodies under rigid bronchoscopy.

6 Indication

[▶ Intended Use, page 18]

7 Contraindication

There are no known contraindications.

8 Special Patient Groups

Does not apply.

9 Product Description

Rotatable, folding bronchoscopy forceps comprising the following parts:

- Handle
- Shaft
- Instrument attachment

Illustration 1: Handle

- A Metal cylinder
- B Thread
- C Twist handle
- D Handle
- E Ratchet (depending on specifications)
- F Release (for handles with a ratchet)
- G Switch-on lever for ratchet (depending on specifications)


Illustration 2: Shaft

- A Swivel nut
- B Ball


Illustration 3: Instrument attachment

- A Jaw
- B Thread


Illustration 4: Cleaning adapter

- A Thread: Connection to the shaft
- B Luer-Lock: Connection to WD

10 Material

- Stainless steel
- Plastic

Not made with natural rubber (latex).

No products made with natural rubber (latex) are used in the production process.

11 Lifetime

[▶ Reprocessing Restrictions, page 20]

12 Shelf Life and Storage

Before first processing:

Store the product in original packaging.

13 Possible Complications and Side Effects

As with underlying disease and procedure.

14 Combining with Other Procedures

Does not apply.

15 Reprocessing

WARNING

- The product is not sterile. Process the product before first and any further application.

This is the only way to ensure the product is germ-free and functional.

Only use the disinfection and sterilization procedures validated for the device and product.

Ensure regular inspection and maintenance of the devices in use.

Adhere to the validated parameters in every cycle.

Observe the local regulations and hygienic rules of the medical office or hospital.

15.1 Reprocessing Restrictions

The product is suitable for multiple reprocessing.

The following factors are relevant for determining the product's lifetime:

- Reprocessing procedure
- Wear and tear
- Chemicals used

Do not use damaged / worn products.

15.2 Cleaning and Disinfection Agents

Do not use cleaning and disinfection agents that contain the following substances:

- Organic, mineral or oxidising acids
- Strong lyes
- Organic solvents (e.g., alcohols, ethers, ketones, benzines)
- Halogens (chlorine, iodine, bromine)
- Aromatic / halogenised hydrocarbons
- Heavy metal salts
- Aldehydes
- Scratching substances

The manufacturer recommends using cleaning and disinfection agents with the following properties:

- Enzymatic cleaning agents to remove protein residue
- Neutral pH value (7.0)

Only use disinfection agents with proven effectiveness.

Only use agents that are compatible with each other and with the devices in use.

Only use agents suitable for cleaning / disinfection of plastic and metals.

Observe information provided by the manufacturer of the cleaning or disinfection agent with regard to concentration and reaction time.

Do not use combined cleaning and disinfection agents.

15.3 Cleaning Preparation

Start: Immediately after product use.

ATTENTION: Remove any foreign bodies from the product.

1. Disassemble the product
[▶ Assemble and disassemble the product, page 23]
2. Clean the product using warm water and standard commercial instruments and cleaning agents.
3. Clean the shaft with an appropriate brush.
4. Remove protein residue using enzymatic cleaning agents.

15.4 Automated Cleaning and Disinfection

Use a tested program for thermal disinfection (max. 95 °C).

Use an appropriate program for the product with a sufficient number of rinsing cycles.

Only use water low in endotoxin (max. 0.25 endotoxin units/ml; e.g., purified water).

Only use sterile or low-germ water (max. 10 germs/ml).

Observe information provided by the WD manufacturer.

The proof of effectiveness of the procedure described here was provided by an accredited testing laboratory using the following materials:

- WD: Getinge Model 88-5
- Disinfection agent: Neodisher medizym (Dr. Weigert GmbH & Co. KG, Hamburg)

When validating your own procedure, please note:

Chemical disinfection may result in disinfection agent residues on the product.

1. Place the product in the WD. The product parts should not come in contact with each other or with other products.
2. Ensure that the jaw of the instrument attachment is open.
3. Attach all lumina of the instrument by means of a suitable rinse adaptor to the rinsing connection of the WD.
4. Start the programme.
5. After successful programme completion: Remove the product.

15.5 Drying

Dry the product immediately after disinfection.

Blow out the tubes.

Only use filtered air for drying.

15.6 Inspection and Function Testing

1. Inspect the product for dirt. Clean and disinfect contaminated products again.
2. Inspect the products for damage (e.g., damaged surfaces, deformations, other mechanical damage). Separate damaged products.
3. Check movable components for ease of movement. Separate the components that do not move with ease.

15.7 Care

1. For maintenance of movable parts use a physiologically safe oil, e.g. paraffin oil in accordance with the European Pharmacopoeia (PH.EUR.) / United States Pharmacopoeia (USP).

15.8 Packaging

Packaging requirements:

- Adequate protection of the product and the sterilization packaging from mechanical damage
- Adequate packaging size (tension-free sealing)
- Sterilization containers: Regular maintenance according to the manufacturer's specifications
- Suitable for steam sterilization (temperature resistant up to at least 137 °C (278 °F), adequate steam permeability)

15.9 Sterilisation

NOTICE

- The temperature in the autoclave must not exceed 137 °C. Otherwise, non-metal components of the product could be damaged.

Procedure	Temperature	Duration (minimum)	Pressure
Fractionated vacuum procedure	134 °C (273 °F)	5 min	2,2 bar

Table 7: Sterilization Parameters

Observe information provided by the autoclave manufacturer.

Ensure that the product is cooled to room temperature before use.

15.10 Storage

After sterilization, store the product dry and dust-free in the sterilization packaging.

Ideally store the products individually in their packaging or in a protective container with compartments.

If the products are stored in drawers: Protect tips with cloth, gauze or small tubes.

16 Application Instructions

WARNING

- Only use intact products from intact sterile packaging or sterilization containers. Only use products that have not expired in terms of their sterilization shelf life.

This is the only way to ensure the product is germ-free and functional.

- Always use the product under visual monitoring.
Otherwise there are risks to the health of your patient.

Ensure the presence of hygienic / sterile conditions needed for the intervention.

Assemble the product only right before use.

Before using the product, carry out a visual inspection and a function test. To do this, examine the product for fractures and cracks. In particular, examine areas such as blades, tips, ends, locks and notches as well as all movable parts. Ensure that the movable parts can be moved.

Disassemble and process the product after use.

16.1 Assemble and disassemble the product

16.1.1 Assemble the product

1. Introduce the instrument attachment into the shaft.

2. Screw on the instrument attachment in the shaft.

3. Ensure that the jaw of the instrument attachment is closed. Pull back the swivel nut.

4. Open the handle completely.

5. Thread the ball at the end of the shaft from the top into the slit in the metal cylinder of the handle.

6. Close the handle. The shaft will be pulled up to the thread.

7. Tighten the swivel nut. **ATTENTION:** Screw the swivel nut firmly against the twist handle.

⇒ The bronchoscopy forceps are now ready for use.

16.1.2 Disassemble the product

To disassemble, detach all parts in reverse sequence.

16.2 Open and close the jaw

1. Open the handle to open the jaw. Close the handle to close the jaw. **ATTENTION:** The opening and closing of the handle does not move the jaw distally or proximally.

2. For handles with a shiftable ratchet: press on the release mechanism to open the handle.

3. For handles with a shiftable ratchet: Move the switch-on lever to position A to activate the ratchet. Move the switch-on lever to position B to deactivate the ratchet.

16.3 Rotate the instrument attachment

1. Rotate the twist handle to rotate the instrument attachment in the desired direction.

17 Instructing the Patient

Does not apply.

18 Aftercare

Does not apply.

19 Maintenance

Check the product during reprocessing for integrity and perform maintenance.

[▶ Inspection and Function Testing, page 22]

[▶ Care, page 22]

20 Disposal

Disposal must be in accordance with national disposal regulations and pursuant to the corresponding risk class.

21 Warranty

The reliability of the product's material and design at the time of shipment is guaranteed. The manufacturer does not know either the diagnosis of the patient or the nature of the application and has no influence on the conditions under which the product is used. The storage conditions after delivery of the product are also beyond the manufacturer's area of responsibility.

Due to biological and individual differences, no product is 100% effective under all circumstances.

Therefore, the manufacturer cannot guarantee a positive effect or the absence of negative effects for product application. The medical staff must use the product on the basis of their medical training and experience, and they are responsible for correct application.

The warranty (repair or replacement) applies only if the product is used in accordance with these Instructions for Use (for instruments, particularly with regard to handling, cleaning, sterilization and maintenance); the warranty period starts on the delivery date.

If you have reason to believe that a new product is faulty, please contact the Customer Service in writing immediately and provide as detailed a description as possible of the fault, the REF (product code), and the LOT (batch code) and/or series number. All allegedly defective products must be returned to us for inspection. Instruments have to be completely cleaned and sterilized, appropriate documentation must be enclosed with the return.

If the manufacturer finds that despite all due care the product was defective at the time of delivery, he will repair the product or replace it promptly. If repair or replacement of the product is not possible, the buyer has the right to cancel the purchase or to reduce the payment, but by a maximum of the purchase price amount.

Additional claims or those not mentioned here due to defect, and other claims regardless of the legal reason, including those based on illegal acts and for compensation of immaterial damages against the manufacturer, his agents, dealers and suppliers, are excluded unless existing law is contrary to the liability exclusion, e.g. in cases of intent or gross negligence or in the event of physical injury.

All claims based on the consequences of non-compliance with the Instructions for Use, including specified indications, contraindications, warnings, instructions, application, storage and off-label use, as well as the consequences of a combination with third-party products are excluded.

Furthermore, all claims that result from the use of products that have expired, or were used despite the obvious damage to the packaging, or resterilized and/or recycled contrary to the Instructions for Use, are excluded.

No one is allowed to change the above conditions, make further warranty or liability declarations, or guarantee any properties that surpass those specified in the Instructions.

The General Terms and Conditions of the manufacturer, which can be accessed at <http://www.novatech.fr/gtc> apply in all remaining instances.

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