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1 Safety Information Marking

WARNING

Non-observance may result in serious injuries, serious deterioration of the general condition or the death of the patient, user, or a third party.

NOTICE

Product damage or other damage may occur in case of non-compliance.

2 Important Safety Information

WARNING

- Before using the product, read the Instructions for Use. Adhere to and save the Instructions for Use.

Otherwise there are risks to the health of your patient.

- Application only by a physician trained in the procedure.

Otherwise there are risks to the health of your patient.

- Do not disassemble or modify the product.

Otherwise there are risks to the health of your patient.

3 Product Codes / REF

[► Product Description, page 12]

4 Package Contents

- 50 x SURVAC™
- 1 x Instructions for Use

5 Intended Use

Removal of foreign bodies, liquids and tissue from the trachea and bronchi.

For use with rigid bronchoscopy.

6 Indication

[► Intended Use, page 12]

7 Contraindication

Does not apply.

8 Special Patient Groups

Does not apply.

9 Product Description

REF	Product name	inner Ø x outer Ø [mm]	Working length [cm]	Material (shaft)
77622	SURVAC™ PTFE 2 x 3 L 60	2 x 3	60	PTFE ¹
77632	SURVAC™ PC 2 x 2,7 L 60	2 x 2,7	60	PC ²
77642	SURVAC™ PC 3 x 4 L 60	3 x 4	60	PC ²
77642SC	SURVAC™ PC 3 x 4 L 30	3 x 4	30	PC ²

1 PTFE = polytetrafluoroethylene, flexible

2 PC = polycarbonate, rigid

- Funnel connector
- Atraumatic ends

Suitable for lateral and axial aspiration (depending on the bronchoscope used).

10 Material

- Connector: Polyvinyl chloride (PVC)
- Shaft: [▶ Product Description, page 12]

Not made with natural rubber (latex).

No products made with natural rubber (latex) are used in the production process.

11 Lifetime

Does not apply.

Single use product.

12 Shelf Life and Storage

For date of expiry, see the product label.

Store the product in unopened original packaging.

13 Possible Complications and Side Effects

There are no known complications or side effects.

14 Combining with Other Procedures

WARNING

- Laser therapy, argon plasma therapy, high-frequency surgery, and other procedures, the effect of which is due to heat: Do not use those methods directly on the product.

Otherwise, injury to the tissue and product damage are possible.

15 Reprocessing

WARNING

- Single use product: Do not reprocess (e.g., clean, disinfect, sterilize), resterilize or reuse the product.

This is the only way to ensure the product is germ-free and functional. Due to the mechanical properties of the product, reprocessing or resterilization could lead to material degradation.

16 Application Instructions

WARNING

- Do not use the product if the packaging or the product is damaged or expired.

This is the only way to ensure the product is germ-free and functional.

Ensure the presence of hygienic / sterile conditions needed for the intervention.

1. Attach the funnel connector of the suction cannula to the suction unit.
2. For lateral suction: Introduce the suction catheter into the lateral port of the bronchoscope.
3. For axial suction: Introduce the suction catheter into the axial port of the bronchoscope.

If suction is weak: Check connection.

17 Instructing the Patient

Does not apply.

18 Aftercare

Does not apply.

19 Maintenance

Does not apply.

20 Disposal

Disposal must be in accordance with national disposal regulations and pursuant to the corresponding risk class.

21 Warranty

The reliability of the product's material and design at the time of shipment is guaranteed. The manufacturer does not know either the diagnosis of the patient or the nature of the application and has no influence on the conditions under which the product is used. The storage conditions after delivery of the product are also beyond the manufacturer's area of responsibility.

Due to biological and individual differences, no product is 100% effective under all circumstances.

Therefore, the manufacturer cannot guarantee a positive effect or the absence of negative effects for product application. The medical staff must use the product on the basis of their medical training and experience, and they are responsible for correct application.

The warranty (repair or replacement) applies only if the product is used in accordance with these Instructions for Use (for instruments, particularly with regard to handling, cleaning, sterilization and maintenance); the warranty period starts on the delivery date.

If you have reason to believe that a new product is faulty, please contact the Customer Service in writing immediately and provide as detailed a description as possible of the fault, the REF (product code), and the LOT (batch code) and/or series number. All allegedly defective products must be returned to us for inspection. Instruments have to be completely cleaned and sterilized, appropriate documentation must be enclosed with the return.

If the manufacturer finds that despite all due care the product was defective at the time of delivery, he will repair the product or replace it promptly. If repair or replacement of the product is not possible, the buyer has the right to cancel the purchase or to reduce the payment, but by a maximum of the purchase price amount.

Additional claims or those not mentioned here due to defect, and other claims regardless of the legal reason, including those based on illegal acts and for compensation of immaterial damages against the manufacturer, his agents, dealers and suppliers, are excluded unless existing law is contrary to the liability exclusion, e.g. in cases of intent or gross negligence or in the event of physical injury.

All claims based on the consequences of non-compliance with the Instructions for Use, including specified indications, contraindications, warnings, instructions, application, storage and off-label use, as well as the consequences of a combination with third-party products are excluded.

Furthermore, all claims that result from the use of products that have expired, or were used despite the obvious damage to the packaging, or resterilized and/or recycled contrary to the Instructions for Use, are excluded.

No one is allowed to change the above conditions, make further warranty or liability declarations, or guarantee any properties that surpass those specified in the Instructions.

The General Terms and Conditions of the manufacturer, which can be accessed at <http://www.novatech.fr/gtc> apply in all remaining instances.

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