
Mode d'emploi
Instructions for Use
Gebrauchsanweisung
Modo de empleo
Istruzioni d'uso
Gebruiksaanwijzing

STERITALC®

F2 / F4 - Flacon - Vial - Flasche - Frasco - Flacone - Flacon

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PF3 - avec kit de poudrage - with poudrage kit - mit Zerstäuberkit - con kit de insuflación - con kit di polverizzazione - met verstuiverkit

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STERITALC® – Medical Talc from Luzenac

REF	STERITALC® Description	Quantity of medical talc	Items/box
16903	F2: Vial, 50 ml	2 g	4
16913	F4: Vial, 50 ml	4 g	4
16833	Spray, with cannula 440 mm	3 g	2
16863	PF3: Vial, 10 ml with balloon and cannula 420 mm (Poudrage Kit)	3 g	2 kits
16983	<i>Accessory for 16863:</i> Vial, 10 ml	3 g	4

The Talc

STERITALC® is non-soluble and induces permanent pleurodesis. It is guaranteed asbestos and latex free. STERITALC® is supplied sterile and has the perfectly adapted and controlled granulometry to minimise the risk of migration across the parietal pleura. STERITALC® has side effects similar to cyclines. Any pain occurring can be reduced by using xylocaine 1% administered directly in the pleura or via the slurry. STERITALC® complies with the European Pharmacopeia.

Indications

Chronic pleurisy, principally malignant.
Spontaneous pneumothorax.
Other cases in which a pleurodesis is indicated.

Contraindications

STERITALC® should not be used if the patient cannot undergo thoracoscopy. In order to avoid systemic migration of the talc, STERITALC® should not be applied after a mechanical abrasion of the pleura. STERITALC® should not be applied in case of pregnancy or breast-feeding.

Dosage

Dosage according to doctor's discretion and in consideration of diagnosis, patient's age and state of health. The dose range varies from 3 to 6 g for pleural effusion and generally does not exceed 2 g for spontaneous pneumothorax. If fluid accumulates in the pleural cavity, such fluid should be drained before talc application.

Precautionary measure

As a precaution, yet another unit of the product should be on standby before the intervention.

Warning (REF 16833)

Do not continuously spray STERITALC® at the same point as there is a risk of congelatio (frostbite).

Instructions for use

- Contents are sterile unless package is damaged or opened.
- Remove contents from inner packaging only immediately before use.
- The respective asepsis standards have to be adhered to when removing the content from the inner packaging.


For the application of STERITALC®, poudrage via thoracoscopy is the method of choice, as it allows for an even distribution onto the parietal and the visceral pleura. Alternatively, STERITALC® (F2 / F4) can be mixed with sterile isotonic solution (NaCl) (> slurry) and applied via a drain or syringe through a chest tube.

Talc induces an irritation resulting in a firm connection between the two pleurae (permanent pleurodesis). An unspecific immune reaction is regarded as normal.

Poudrage

Spread the talc homogeneously inside the pleural cavity. **The aim is to produce a “mist” inside the pleural cavity.**

- **STERITALC® F2 / F4:** Connect the vial to a nebulizer or decant the contents of the vial into the nebulizer.
- **STERITALC® Spray:** Shake the spray can. After removal of the cap, fasten the actuator with the cannula securely onto the top of the spray can. Before application on the patient, depress the actuator button once shortly to make sure the actuator being fixed correctly and the spray working properly. Hold the spray can in an upright position. **Spray evenly and intermittently to the areas indicated for talc application. Do not continuously spray at the same point.**
- **STERITALC® PF3:** Fully remove the aluminium tear-off cap from the STERITALC® vial. Remove the plug from the vial. Make sure that the cannula is attached to the white coupling. Firmly press the coupling onto the vial until the coupling “clicks” onto the vial top. This guarantees a tight fit and seal necessary for full functionality of the device. Attach the balloon to the coupling with the luer connector. Make sure all components are correctly assembled.

About 1 to 2 minutes are necessary to apply 3 g of talc, pressing the balloon regularly with a maximum frequency of 1 injection per second. The manufacturer has made allowance for a small quantity of talc remaining inside the vial after the procedure, which is normal and does not affect the efficiency of the talc application.

Introduce the cannula into the trocar and start pulverization (F2 / F4 / PF3 / Spray). The first impulses should be gentle so as to not apply too large a quantity of talc. In order to achieve a homogeneous distribution of the talc in the pleural cavity, repeat procedure changing the orientation of the cannula and making sure that the talc does not adhere on the inner walls of the vial / can, which is achieved by a slight swinging motion of the vial / can.

Warning: The cannulas provided with STERITALC® Spray and STERITALC® PF3 have been specially designed to minimize the risk of a lung perforation and to optimize talc flow. If necessary, the cannula can be shortened by cutting it precisely and free of burrs with a new scalpel. The use of any other cannula system is the physician's responsibility as it can affect the quality and duration of talc application and may cause injury to the patient.

With the application of talc, air is introduced into the pleural cavity. Therefore the trocar used for application must not be closed tightly. Follow instructions for use of the trocar.

Slurry

- **STERITALC® F2 / F4:** STERITALC® can be used in a slurry form. Novatech recommends the use of 10 ml sterile isotonic solution (NaCl) for each 1 g of talc. Pain can be reduced by adding 1% xylocaine. A slurry must not be prepared in advance. Use mixture immediately. Shake mixture thoroughly to avoid agglutination! Apply via a drain or a 50 ml lavage syringe and inject in the pleural cavity through a chest tube. Changing the patient's position contributes to a homogeneous distribution of the suspension and is recommended.

Patient Information

Prior to x-ray examination the radiologist must be informed of the talc pleurodesis. As the talc remains inside the pleural cavity it could be wrongly identified as a nodule, tumor or cyst.

Literature

Please also refer to the literature published about the use of STERITALC® listed at <http://www.novatech.fr/pneumologie/steritalcr.html>

Additional Information

Disposal according to national regulations for waste and / or - after patient contact or contact with contaminated material - to regulations for contaminated waste.

Products and packaging are free of latex.

STERITALC® Spray:

Can is pressurized. Do not puncture or burn.

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