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
### Intended use

The TONN placement kit is designed to fold and insert Dumon® silicone stents without kinking into bronchoscopes and to place them in cases of benign and malignant tracheobronchial stenoses.

### Instructions for use

Place the basic block onto the folding system and lubricate generously (1). We recommend using a lubricant developed by NOVATECH (pharmaceutical grade). **Lubrication is essential in order to facilitate placement and to avoid stent damage.** Lubricate the stent, too (2).

Position the stent centrally on the folding system (3). Insert the clamping rod by inserting its end into the provided hole in the basic block (on the side with the NOVATECH logo). Fold the stent by pushing down the clamping rod (4). Block the clamping rod in the basic block by rotating the rod one quarter.

Close the folding system ensuring that the stent is pushed inside completely (5). Lock the folding system, pushing the slider with the coloured button

first to the right side, then to the left. The folding system must be closed completely (6).

Remove the clamping rod together with the folding system by rotating the rod one quarter. Pull the clamping rod out of the folding system (8). Note: In case of a lack of lubrication the stent might come out along with the clamping rod. In this case, simply start the procedure from the beginning.

Lubricate the inside of the introducer tube generously (9). Insert the introducer tube into the folding system (see mark on the folding system, 10).

By assuring that the NOVATECH logo is on top of the introducer tube (11), you can identify the position of the stent at any time. This is especially beneficial when Y-stents are placed (see „Placement of Y-stents“). Insert the loading rod into the folding system (see mark on the folding system, 12). Hold the folding system with the inserted introducer tube and loading rod in a vertical position.

Apply light pressure to the loading rod, carefully pushing the stent inside the tube (13). The loading rod cannot be pushed beyond the security point (14). Withdraw the introducer tube from the folding system. The stent is positioned inside the introducer tube (15).

Introduce the pusher into the proximal end of the introducer tube. The graduation on the pusher indicates the exact length of the stent inside

the introducer tube (16). By applying pressure onto the pusher, the stent is easily released (17, 18) through a STORZ bronchoscope.

## Placement of Y-Stents

The position of the main branches can be checked before placement and the risk of damaging such a stent is minimised. Lubricate the applicator system as described above. Place the Y-stent in the folding system with its right branch on the upper side of the folding system (19). Fold the stent with the help of the clamping rod (20) as described above. Close the folding system locking the stent inside (21). Important: Ensure that the NOVATECH logo is on top when inserting the introducer tube (22).

This ensures that the right branch of the stent will be placed on the right hand side of the logo and will therefore be inserted in the main right bronchus (23, 24). Keeping the NOVATECH logo on top assures correct placement of the stent within the trachea before placing it in the main carina. This result is obtained with all Y-stents and tracheal stents with diameters of 18 and 20 as those stents have no studs on the posterior part. The red TONN/NOVATECH system enables the user to place large Y-stents.

## Warning and precautionary measures

Use only cleaned and sterilised instruments from intact sterile packaging or sterile containers.

Effective cleaning and disinfection is a vital prerequisite for effective sterilisation.

Within the scope of your responsibility for the sterility of the instrument during its application, please take the following into account:

- Only apply equipment- and product-specific cleaning / disinfection and sterilisation procedures that have been adequately validated.
- Regularly service and check the equipment used (cleaning and disinfection device, sterilizer) and
- ensure that the validated parameters are maintained during each cycle.

Please only use defect-free products as a matter of principle and give each instrument a brief visual and functional check before using it on the patient. Do not use in case of doubt and contact your NOVATECH service partner.

Please also note the statutory and hygiene provisions of the doctor's surgery / hospital applicable in your country, particularly with regard to the different guidelines on the effective inactivation of prions.

## Reprocessing restrictions

The TONN placement kit and STORZ bronchoscopes are suitable for multiple reprocessing in accordance with the instructions set out below. The end of product life is normally determined by wear and tear and damage through use. Damaged and worn instruments or instrument parts have to be sorted out and discarded. The responsibility for the further use of damaged and / or contaminated devices lies solely with the user.

When selecting the cleaning or disinfecting agents, please make sure that they do not contain the following ingredients:

- organic, mineral or oxidising acids (corrosive liquids); a neutral / enzymatic cleaner is recommended
- strong lyes

- halogens (chlorine, iodine, bromine)
- organic solvents (e.g. alcohols, ethers, ketones, benzines)
- aromatic / halogenated hydrocarbons
- heavy metal salts.

Never clean the instruments, sterilisation trays and sterilisation containers with metal brushes or wire wool.

Instruments, sterilisation trays and sterilisation containers may only be exposed to a maximum temperature of 141 °C (286 °F) unless stated otherwise by the manufacturer.

### Place of use

Please ensure that contaminated instruments and products are collected separately and not returned to the instrument trays during the application to avoid further contamination of the loaded instrument trays.

### Storage and transport

No special requirements apply for storage and transport.

The instructions for the cleaning of the instrument set out below have to be performed within 2 hours after the use of the instrument.

### Cleaning preparations

For cleaning, instruments must be disassembled first:

Open the TONN placement kit and ensure that the folding system is removed from the basis block during cleaning / sterilisation.

Remove all accessories from the STORZ bronchoscope and clean the disassembled components separately.

Pre-treatment measures must be performed before automatic and manual cleaning and disinfection. Please use running water or a disinfectant to remove gross contamination. The applied disinfectant should be aldehyde-free, its effectiveness should be documented and it should be suitable for and compatible with the instruments. To manually remove contaminations, please always use a soft brush or a clean soft cloth which you only use for this purpose. Never use metal brushes or wire wool!

Rinse all lumina of the instrument at least five times - more often if required - by means of a single-use syringe (syringe volume at least 10 ml). In addition, all lumina of the bronchoscope should be carefully brushed.

Please note that the disinfectant used for pre-treatment only serves to protect the person in charge of the processing. It cannot replace subsequent disinfection procedures after cleansing.

The bone dust collector is reassembled under aseptic conditions directly before the intervention. In this connection, please note the instructions under "Assembly before use".

### Automatic cleaning process

Please note the following when selecting the cleaning and disinfection device:

- The effectiveness of the cleaning and disinfection device must always be documented (in this respect, please comply with the statutory provisions of your country).
- A certified thermal disinfection program must be used (for a minimum of 10 min at 93 °C).
- It must be ensured that the applied program is suitable for the instruments and contains an adequate number of rinsing cycles.
- It must be ensured that only sterile water or water with low germ content (max. 10 germs/ml) and water with low endotoxin content (e.g. purified water / highly purified water with max. 0.25 endotoxin units/ml) is used for mid-cycle or final rinses.
- It must be ensured that only filtered air is used for the drying process and
- that the cleaning or disinfection device is regularly serviced and checked.

When selecting the cleaning agent system, it must be ensured that it is suitable for the cleaning of instruments made of metal and plastic materials, and compatible with the cleaning and disinfection device that you use.

The instructions issued by the cleaner / disinfectant manufacturer or the cleaning and disinfection device manufacturer with regard to concentrations and the use of the products must be adhered to at all times.

Procedure:

1. Place the already disassembled instrument into the cleaning and disinfection device. Please ensure that the instrument parts do not come into contact with each other. Attach all lumina of the instrument by means of a suitable rinse adaptor to the rinsing connection of the cleaning and disinfecting device.
2. Start the program (disinfection for a minimum of 10 minutes at 93°C).
3. Remove the instrument parts from the cleaning and disinfection device after the successful completion of the program.

Please execute the following steps immediately if possible:

4. Dry the instrument parts if required.
5. Check the instrument parts for any damage.
6. If required, sort the cleaned and disinfected instruments in suitable the instrument trays.
7. Package the instrument parts / the loaded instrument trays.

Regardless of the above indicated restrictions, please find a list of agents suitable for automatic cleaning and disinfection on page 23.

In case you wish to validate your own process, please consider the risk of disinfectant residues on the product during chemical disinfection. In case of doubt, please contact your NOVATECH service partner.

### Manual cleaning process

An automatic cleaning process (cleaning and disinfection device) as described above should, if possible, be applied for cleaning and disinfection. Due to its significantly reduced effectiveness and repeatability, a manual process - including the use of an ultrasound bath - should only be applied if an automatic cleaning method is not available.

When selecting the cleaning and disinfecting agents, please ensure the following:

- The cleaning and disinfecting agents should be suitable for use on instruments made of metal and plastic materials.
- The cleaning agent - if applicable - should be suitable for ultrasonic cleaning (no foam formation) and
- a disinfecting agent with a certified effectiveness and compatible with the applied cleaning agent should be used (in this connection, please also note the statutory provisions in your country).

Combined cleaning / disinfecting agents should not be used if possible.

The concentrations and residence times stated by the manufacturer must be maintained without fail. Please only use freshly prepared solutions, sterile water or water with low germ content (max. 10 germs/ml) and water (e.g. purified water / highly purified water) with low endotoxin content (max. 0.25 endotoxin units/ml). Use only filtered air for the drying process.

Procedure:

1. Place the already disassembled instrument into the cleaning bath for the stipulated residence time so that all parts are sufficiently covered (possibly with ultrasonic support or careful brushing using a soft brush). Please ensure that the instrument parts do not come into contact with each other. Rinse all lumina five times at the beginning and five times at the end of the residence time by means of a single-use syringe (syringe volume at least 10 ml) and a suitable rinsing adapter.
2. Subsequently remove all instrument parts from the cleaning bath and thoroughly rinse them with water at least three times. Rinse all lumina of the instruments at least five times by means of a single-use syringe (syringe volume at least 10 ml) and a suitable rinsing adapter.

3. For disinfection, place the disassembled, cleaned and inspected instrument parts into the disinfecting bath for the stipulated residence time so that all instrument parts are sufficiently covered. Please ensure that the instrument parts do not come into contact with each other. Rinse all lumina of the instrument by means of a single-use syringe (syringe volume at least 10 ml) and a suitable rinsing adapter five times at the beginning and five times at the end of the residence time.
  4. Subsequently remove the instrument parts from the disinfecting bath and thoroughly rinse them with water at least three times. Rinse all lumina of the instrument at least five times by means of a single-use syringe (syringe volume at least 10 ml) and a suitable rinsing adapter.
- Please execute the following steps immediately if possible:
5. Dry the instrument parts if required.
  6. Check the instrument parts for any damage.
  7. If required, sort the cleaned and disinfected instruments in suitable the instrument trays.
  8. Package the instrument parts / the loaded instrument trays.

Regardless of the above indicated restrictions, please find a list of agents suitable for manual cleaning and disinfection on page 24.

### Disinfection

The applied disinfection process is described in more detail in the “Cleaning” section.

### Drying

If further drying is required, please ensure that it is performed in a clean location to avoid recontamination of the already cleaned and disinfected parts.

### Maintenance

Instrument oils should, if possible, not be used. If their use is required, then it should be ensured that only instrument oils (pharmaceutical white oil / surgical lubricants) are used which are licensed for use in vapour sterilisation and which have verified biocompatibility, taking into account the maximum sterilisation temperature. Use the instrument oil only on mechanically loaded and movable parts (large-scale use to prevent corrosion is not necessary and not permitted). We recommend using a lubricant developed by NOVATECH (pharmaceutical grade).

### Inspection and functional testing

Please perform a visual check for corrosion, damaged surfaces, deformations, chipping and contaminations following the disinfection process. Eliminate damaged and not fully functioning instruments or instrument parts and repeat the cleaning and disinfection process in the case of contaminated instruments. Check movable components for their functional capability and free movement.

If required, get in touch with your NOVATECH service partner.

### Packaging

Package the instruments or the sterilisation trays with the instruments using single-use sterile packaging (single or double packaging) and / or sterilisation containers. To do this the following requirements must be met:

- They conform to the statutory provisions of your country.
- They are suitable for vapour sterilisation (with a temperature resistance of at least 141°C (286°F) and adequate vapour permeability).
- They adequately protect the instruments / sterilisation packages against mechanical damage.
- They are adequately sized to ensure that the sealed seam is not under mechanical strain.
- They are regularly maintained in accordance with manufacturer’s instructions (sterilisation containers).

## Sterilisation

Vapour sterilisation must be applied to sterilise the TONN placement kits and the STORZ bronchoscopes. A fractionated vacuum procedure or gravitation procedure (with adequate product drying) is to be applied.

Please ensure that the sterilizer conforms with the statutory provisions of your country. A sterilisation temperature of 141 °C (286 °F) may not be exceeded. A sterilisation time (exposure time at sterilisation temperature) of at least 5 min<sup>1</sup> at 134 °C (273 °F) must be adhered to.

Hot-air sterilisation, X-ray sterilisation and sterilisation with formaldehyde or ethylene oxide have so far not been validated. If required, please contact your NOVATECH service partner for further information.

## Storage

Ensure that following the sterilisation process the instrument parts in their sterilisation packages are stored dry and dust-free. The replacement filter media BM90KNMS12 must be stored at room temperature in a dark place.

## Assembly before use

The TONN placement kits and the STORZ bronchoscopes are reassembled under aseptic conditions directly before the intervention.

## Additional information

The instructions stated above were validated by NOVATECH as being SUITABLE for the preparation of the processing of a medical device before application. The person in charge of the processing of the product is responsible for ensuring that the actual process leads to the desired results using the equipment, materials and staff in the processing unit. This normally requires process validation and routine monitoring. In addition, any deviation from the instructions provided by the person in charge must be evaluated carefully for effectiveness and possible detrimental effects.

## Contact to manufacturer

Please contact your local service partner or get in touch with us via [www.novatech.fr](http://www.novatech.fr) if you have any queries with regard to the product, any feedback or suggestions for improvement.

## Return consignment

N.B.: An instrument must be cleaned, decontaminated and sterilised before it is sent back to us for processing / repair.

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