

# FleXNeedle<sup>®</sup>

## Multipurpose 18g Coring Needle



The logo for BRONCUS features a stylized yellow and blue circular graphic to the left of the word "BRONCUS" in a bold, blue, sans-serif font.

BRONCUS<sup>®</sup>

The logo for INTERVENE MEDICAL features a stylized blue and grey graphic to the left of the word "INTERVENE" in a bold, blue, sans-serif font, with "MEDICAL" in a smaller, grey, sans-serif font below it.

INTERVENE  
— MEDICAL —

The **FleXNeedle** is designed to enhance accessibility in hard to reach lesions in the peripheral segments of the lung.



#### ENHANCED FLEXIBILITY

- ~ 70% shorter hypo tube length provides enhanced flexibility & accessibility to tortuous airway anatomy

#### GREATER CONTROL

- Ergonomic, single-handed needle actuation
- Standard Luer Fitting: Compatible with wide range of aspiration syringes & stopcocks
- Adjustable depth gauge provides penetration of 5 to 20mm for precision and control

#### OPTIMIZED BIOPSY

- Proprietary coring tip cuts and captures core specimens for histology and cytology
- Needle internal lumen is 0.95 mm enabling a large sample collection

## Ordering information

Catalog Number	Needle Gauge	Penetration Depth (mm)	Maximum Catheter OD (mm)	Catheter Length (mm)
10005-2	18	0-20	1.9	1400

**With Access Comes Possibility™**

### Broncus Medical, Inc.

125 Nicholson Lane, San Jose, CA 95134

Contact: For additional information please contact Broncus at 1-877-428-1600, 1-650-428-1600

or sales@broncus.com

[www.broncus.com](http://www.broncus.com)



Intervene Medical Pty Ltd - Improving Patient Outcomes  
 Head Office: U1/22 Payneham Rd, Stepney, SA 5069, AU  
 Phone: 08 6444 9949  
 Email: [customer.service@intervenemedical.com.au](mailto:customer.service@intervenemedical.com.au)